




The Rightway<sup>®</sup> System

Ultracare featuring The Rightway<sup>®</sup> System offers a superior choice for the prevention of DVT.



Two-Way Stretch Fabric

Alignment Indicator



Rightway Ankle Crest

Prevents skin shearing, irritation and excess material at top of ankle



Rightway Heel Guides

Assists with proper placement on foot



Rightway Heel Pouch

Oversized for added comfort

Top View Inspection Port



## UltraCARE® Anti-Embolism Stockings Features and Benefits

### Patient Comfort and Safety

**Anti-Microbial Finish** Aids in odor management.

**Ribbed Ankle Crest** Prevents excess fabric abrading against the skin.

**Oversized Heel Pouch** More room & improved comfort.

**Comfort Panel** Accommodates changes in limb size and shape.

**Top-View Inspection Port** Better patient comfort when ambulating.

**Rapid Moisture Transfer** Moves fluids and body oils from skin to surface of the fabric quickly to speed drying; launders better.

**Two-Way Stretch Fabric** Allows stocking to conform to shape of limb.

**Premium Elastic Top** Keeps stocking in place without causing tourniquet effect.

**Outward Facing Seams** Reduces skin irritation and/or depressions as limb size changes.

**Latex-Free** Safe for latex-sensitive skin.

### Clinician Focused

**Color-Coded Heel Guides** Assists in proper placement and provides size identification.

**Ribbed Ankle Crest** Prevents excess fabric abrading against the skin.

**Oversized Heel Pouch** Eases donning and doffing while assisting proper placement.

**Comfort Panel** Allows more stretch to position monitors or check incision sites.

**Color-Coded Top-View Inspection Port** Conveniently access toes and top of foot and identify size.

**Alignment Indicator** Visible guide to prevent twisting during application.

**Color-Coded Packaging** High visibility colors for easy size identification and selection.

**20-30mmHg Graduated Compression** helps improve blood flow and aids in clot prevention.

### Available Sizes and Styles

**Over 40 Sizes To Choose From** Knee, thigh and waist length.

For more information, visit [www.DVTSTAR.com](http://www.DVTSTAR.com)



## UltraCARE® Anti-Embolism Stockings

Short	Regular	Long	Style	Size	Color	Qty.
—	853-00	—	Knee	X-Small	Pink	12 PR
851-01	853-01	858-01	Knee	Small	Yellow	12 PR
851-02	853-02	858-02	Knee	Medium	Orange	12 PR
851-03	853-03	858-03	Knee	Large	Blue	12 PR
851-04	853-04	858-04	Knee	X-Large	Green	12 PR
851-05	853-05	858-05	Knee	2X-Large	Grey	12 PR
—	853-06	—	Knee	3X-Large	Purple	12 PR
863-01	868-01	873-01	Thigh	Small	Yellow	12 PR
863-02	868-02	873-02	Thigh	Medium	Orange	12 PR
863-03	868-03	873-03	Thigh	Large	Blue	12 PR
863-04	868-04	873-04	Thigh	X-Large	Green	12 PR
863-05	868-05	873-05	Thigh	2X-Large	Grey	12 PR
—	868-06	—	Thigh	3X-Large	Purple	12 PR
—	83004-01	83007-01	Waist	Small	Yellow	6 PR
—	83004-02	83007-02	Waist	Medium	Orange	6 PR
—	83004-03	83007-03	Waist	Large	Blue	6 PR
—	83004-04	83007-04	Waist	X-Large	Green	6 PR